





Six-Dimensional Activation
6D TAPE®

Small handles – great benefits!



*Six-Dimensional
Activation*





6D Tape Inc.



6D Tape Inc. is a Finnish health technology company founded in 2016. We focus on the development of medical devices and products for physical therapy, lymphatic therapy and active self-care.

Our mission is to produce next-generation solutions for all who need faster healing and an improved quality of life.

Global Headquarter is in Sipoo, Finland.

Current product portfolio

6D Tape products: CE MDR 2017/745/EU: Class 1; USA FDA Class 1, 510(K) Exempt, KGX, #21 CRF 880.5240, Listing #D442239, Patented.

6D Action Device products: EU CE, EN 60335-1:2012 +A11:2014; USA FDA Class 1, 510(K) Exempt, ISA, #21 CRF890.5660, Listing #D442242.

Small handles - great benefits!
6D TAPE®

Small handles - great benefits!



Contacts

Tapani Taskinen

CEO, D.Sc. in Industrial Engineering & Management

Tel. +358 50 544 8444

tapani.taskinen@6dtape.com

Jussi Karkela

Sales & Clinical Director, B.Sc. in Physiotherapy

Tel. +358 44 276 3987

jussi.karkela@6dtape.com

Head Quarter Address

6D Tape Inc.

Iso Kylätie 47

04130 Sipoo, Finland

www.6dtape.com

Part 1: 6D Tape – Background – Theory – Principles



6D Tape: Core features & technologies

6D Tape has many therapeutic uses

6D Tape contraindications

6D Tape treatment activates the lymphatic system

6D Tape feasibility study

6D Tape treatment examples & results



Part 2: How to treat with 6D Tape?

Treatment protocol for 6DTape manipulation

Watershed areas and directions of flow

6D Taping in practice

Treatment videos / training with videos

summary

6D Tape taping Video



Treating the lower back symptoms with 6D Tape:

1. Start: activating the left venous angle
2. Treatment of lower back symptoms
3. End: activating the left venous angle



6D Tape core technologies

1 – 6D manipulation



Six directions are along the X, Y and Z axis, including rotation around the axis.

Controlled, painless and safe manipulation



Six-Dimensional Activation

Six-Dimensional Activation
6D TAPE[®]

6D Tape core technologies

2 – skin-friendly adhesive

The left picture shows how the 6D Tape adhesive is gentle on the skin when the tape is removed. The picture on the right shows how the adhesives of regular kinesiology tapes and other therapeutic tapes damage the protective top layer when removed.

Six-Dimensional Activation
6D TAPE®


Small handles – great benefits!

Stratum Corneum Protection System

6D Adhesive is less damaging to the stratum corneum when it is removed

6D Tape Adhesive

Conventional adhesive




Less stratum corneum is removed.

Skin damage may be caused by removal of stratum corneum

Microscopic view of "adhesive surface after removal"

6D Tape core technologies # 3 – excellent adhesion



6D Tape Adhesive contours to the surface of the skin, as shown in the picture on the left.



Regular tapes use strong adhesives that will not conform to the skin's surface, as seen in the right picture. This is why regular tape will damage the top layer of the skin when removed

6D Tape adhesive technology is skin-friendly



6D Tape core technologies # 4 – skin friendly base material



1. A unique, elastic and skin-friendly kinesiology tape that is 98 percent cotton



6. No latex or rubber is used

2. The 6D Tape is elastic between the handles

3. 6D Tape will also breathe, and it is designed to be worn for 2–10 days.

4. waterproof, making it safe for saunas, showers, swimming and so on.

5. Hypoallergenic and suitable for use on highly sensitive skin..

Six-Dimensional Activation
6D TAPE®

Small handles – great benefits!

6D Tape core technologies # 5 – skin friendly handle material



1. The 6D handles are soft, as they are made with 100 percent polypropylene fibers that are used in baby diapers, for example.

2. A good grip, as the material is not flexible

3. Handle will also breathe, and it is designed to be worn for 2–10 days.

4. waterproof, making it safe for saunas, showers, swimming and so on.

5. Hypoallergenic and suitable for use on highly sensitive skin..

6. No latex or rubber is used




Six-Dimensional Activation
6D TAPE®

Small handles – great benefits!


6D Tape core technologies

6 – active self-care 24/7





Six-Dimensional
Activation



Six-Dimensional Activation
6DTAPE®

Small handles – great benefits!

01/14/2025 © 6D Tape Oy

Indication Segments

Edema, Swelling, Scar, Wounds, Fascia, Fibrosis, Muscles and Pain




Chronic Edema




Pains & Aches

Six-Dimensional Activation
6D TAPE®
Small handles - great benefits!



Pressure Wounds



Pre- and Post-Operative



Acute & Chronical Injury

6D Tape indications



- *Scars and fibrosis*
- *Swelling and pain*
- *Acute and chronic injuries*
- *Back pain, herniated discs*
- *Sprains and injuries to muscles and ligaments*
- *Neck pain*
- *Shoulder injuries and pinched nerves*
- *Joint problems*
- *Plantar fasciitis*
- *Medial tibial stress syndrome, or shin splint*
- *Golf and tennis elbow*
- *Postoperative rehabilitation*
- *Lymphedema*
- *ITB syndrome, hip bursitis*
- *Osteoarthritis*
- *Bruising*
- *Sports injuries*
- *Etc.*



6D Tape absolute and relative contraindications



* Contraindications of physiotherapy ,lymphatic therapy and massage

- Acute infections – no tape on the site
- Fresh blood clots
- Ongoing cancer treatment *
- Severe heart failure
- Do not put 6D Tape on open wounds
- Do not put 6D Tape on stiches
- No tape under garments




Small handles – great benefits!





Lymphatic system






Lymph flow is part of the body's immune system.

Different parts of the body contain lymphatic tissue that protects the body from microorganisms and harmful substances. Lymph fluid transports lymphocytes, or white blood cells, to the lymph nodes, where the lymphocytes destroy harmful substances and eliminate pathogens.

The function of the lymphatic system is to return the substances belonging to the lymphatic load from the extracellular fluid, or tissue fluid, back into the bloodstream. Such substances include water, proteins, dead cells, fatty acids, cell metabolites, bacteria, and viruses.

The key is to understand that if proteins cannot exit the extracellular space, for example as a result of tissue damage, they absorb fluid via osmosis, which results in the swelling of the tissue.

Why 6D Tape works?



Six-Dimensional Activation
6D TAPE®

Small handles - great benefits!

Why 6D Tape works?



Proteins


Close-up of a lymphatic capillary

Intracellular space



➤ **Filamentti proteiini**

Why 6D Tape works?




Lymph node

Lymphatic collector

Proteins

Close-up of a lymphatic collector




Six-Dimensional Activation
6D TAPE®

Small handles – great benefits!



6D Improves Lymphatic Circulation

- ❑ Active manual self-treatment
- ❑ Reduces swelling and pain
- ❑ Releases scars and fibrosis
- ❑ Accelerates healing
- ❑ Tool to support the work of the physiotherapists and other professionals




Scientific Proof with Imaging



Proved effect on
lymph flow:

6D Tape® (right)
vs
MLD (left)



6D Tape feasibility study



The research question was
**‘Does 6D Tape have a local
effect on the flow of lymphatic
fluid in a normal lymphatic
system?’**

The study was coordinated by:

Jean-Paul Belgrado, PhD, Email: belgrado@ulb.ac.be

Lymphology and Rehabilitation Research Unit's

Faculty of Motor Skills Sciences, Laboratory of Sciences of
Motricity

Université Libre de Bruxelles

CP 640

Route de Lennik 808

1070 Brussels

Belgium

Watch a video on the feasibility study on the 6D Tape website:

<https://www.6dtape.com>

6D Tape feasibility study


The **research method** chosen was contrast medium imaging. The medium was injected into soft tissue and its passage and rate of travel were recorded on video in real time. The visualization was created by using two cameras: one recording the skin and the other recording the movement of the medium under the skin.

The upper right image shows the contrast medium being injected into the tissue. The lower picture shows the measuring lines between which the flow of the contrast medium is examined. The interval between the measuring lines was 50 mm



Screenshot of a video taken with an infrared camera:
Fluorescent agent is injected into the tissue


The distance between the lines of measurement of the fluorescent agent test was 50 mm.



A: Rest



C: Manual lymphatic drainage



B: Voluntary mobilization



D: Mobilization with a single 6D Tape handle



E: Mobilization with five handles of the 6D Tape

Near Infrared Fluorescence Lymphatic Imaging (NIRFLI)
of lymph flow

vs

Manual lymphatic drainage MLD

6D Tape with mobilization handles



6D Tape feasibility study

Experimental condition:

- **Condition A:** Rest
- **Condition B:** Voluntary mobilization (flexion – extension)
- **Condition C:** Manual lymphatic drainage (MLD, Fill & Flush Method)
- **Condition D:** Mobilization with a single 6D Tape handle
- **Condition E:** Mobilization with five handles of the 6D Tape



Condition D 6D Tape Feasibility Test



The 6D Tape strip that is manipulated is placed over the right side of the hand.

6D Tape feasibility study

6D Tape feasibility
test video short
version


<https://youtu.be/Xz-O6t7BxwY>

Scientifically Tested with NIRFLI-imaging



6D Tape is the first and only therapy tape that has positive scientifically proven results with Near Infrared Fluorescence Lymphatic Imaging for improving fluid and lymph flow. In the Table below 6D Tape treatment is compared with manual lymphatic drainage, muscle mobilization and rest condition.

	Rest	Muscle Mobilization	MLD, Fill&Flush	6D Tape / One Handle Strip			6D Tape / Five One Handle Strips		
Variable/Condition	A	B	C	D 1	D 2	D 3	E 1	E 2	E 3
Presence of Lymph Bolus Traffic	Very limited	Small sign	Clear sign	Clear sign	Clear sign	Clear Sign	Clear sign	Clear sign	Clear Sign
Increase of Lymph Bolus Traffic	No	Very limited	Pulse in each manipulation	Continuous but low	Continuous but low	Continuous but low	Continuous, rapid in the end	Continuous, rapid in the end	Continuous, rapid in the end
Sequence Length [s]	300 sec	74 sec	42 sec	70 sec	187 sec	148 sec	29 sec	47 sec	32 sec
Manipulation Technique	Rest	Wrist Movement	Fill & Flush	Pull, Rotate, Horizontal Move			Pull, Push, Rotate, Horizontal Move		



Conclusion: 6D Tape treatment promotes fluid flow from the extracellular space into the lymphatic capillaries and further into the superficial lymphatic vessels.

The Study is published in the journal *Lymffis*, issue 2/2019, in Finnish. The study is also available in English and in German.

For more information about the research look at <https://www.6dtape.com/research-results/> and download the research paper.



Small handles – great benefits!


Our Customers Love Us



“This tape has proven useful in self-management to help stimulate lymph drainage more pro-actively.

It is also extremely effective at helping release tight scar tissue – which again helps reduce localized congestion.”

Excerpt from the book Your Lymphoedema – Taking Back Control, Lymphedema Support Network, UK.



Small handles – great benefits!
6DTAPE[®]

Small handles – great benefits!
Six-Dimensional Activation

Read more professional user experience stories at <https://www.6dtape.com/blog/>

Small handles – great benefits!

Benefits of one hour edema treatment with 6D Tape




Before



During

1 hour



After

Small handles - great benefits!
6D TAPE®

Small handles - great benefits!

6D Tape – A Breakthrough in Advanced Wound Care and Scar Management



- In the world of wound care and rehabilitation, the quest for effective and innovative solutions is ongoing.
- We are excited to introduce **6D Tape, a game-changer in managing edema, swelling, scar tissue, and pain.**
- Combined with surgical post-operative film dressings that inhibit bacterial growth (like Leukomed Sorbact), 6D Tape provides a holistic strategy for managing wounds and scars. This is particularly effective for post-C-section recovery as well as orthopedic wounds and scars that arise from both trauma and surgery.



Small handles - great benefits!
6D TAPE®

Small handles - great benefits!

Treatment of a Scar Tissue in the Wrist



One week treatment with 6D Tape

Not treated with 6D Tape



Benefits for the patient – as told by 6D Tape customers

A customer of Physiotherapist and Lymphatic Therapist Christian Samp underwent knee joint surgery in October 2019 in Germany.

The customer now has issues with swelling and scarring as well as a limited range of motion.

The video shows the patient treating himself with 6D Tape – a one-handle piece of tape is placed on the venous angle in the supraclavicular fossa.

The person works with 6D Tape 5 times a day.

The result is a 2 cm smaller knee circumference in one day!


Knee 6D taping will be discussed in more detail in the treatment section of this presentation later.

Benefits for the patient – as told by 6D Tape customers



Said by our customers: “Background: A 30-year-old breast cancer survivor, with breast tissue removed from both sides and replaced with silicone implants.

In this example, a 30-year-old woman shares her experience of 6D Tape treatment after breast cancer surgery. Her breast tissue was replaced on both sides with silicone implants. Subdermal scar tissue has developed in the breast and chest area after the surgery, the skin of the breasts has attached to the implants, and the subdermal tissue and fasciae have constricted considerably. This caused a sensation of pressure and tightness in that region, along with shoulder and neck pain, as well as a restricted range of motion.



Observations and experience of 6D taping:

- **The treatment with 6D Tape provided welcome additional help in treating the issues and made a big impression.**
- **With the tapes, I was able to treat the problem areas more easily than without them.**

Six-Dimensional Activation
6D TAPE®

Small handles – great benefits!

- Even the tighter skin can be treated with the tapes, while without the tapes, you couldn't get a good grip on the skin. In particular, lifting the skin and pulling in different directions has been an added benefit of the tape treatment.
- After treatment, the tissue is softer and more flexible. My skin is becoming more elastic week after week.

Benefits for the patient – as told by 6D Tape customers



- Manual handling of areas takes physical strength and can be difficult, but the handles are easy to hold and they can be used to mobilize the area in multiple directions.
- When the tapes are placed on the area to be treated, they also remind me to repeat the treatment more often.
- Placing the tapes was easy. The tape has not caused any skin problems, not even any irritation or redness, although normally even band aids left in place for a few hours cause irritation for me.
- I have also used 6D Tape for treating sore knees, for example. I've had long-lasting pain there that was reduced after just a week of "tape wiggling" and other stretching.



With 6D Tape, treatment is more efficient and easier compared to treatment without tape. It's also inexpensive compared to physiotherapy visits."

Six-Dimensional Activation
6D TAPE®

Small handles – great benefits!

6D Action Techniques Video



Pulsation

**Lift
and twist**

Gliding



Part 2: How to treat with 6D Tape?

Treatment protocol for 6D manipulation

Watershed areas and directions of flow

6D Taping in practice

Treatment videos / training with videos

summary

6D Tape treatment protocol




1. Getting started

Always start the treatment from the lymph node area located in the venous angle. (one minute, 10-15 manipulations)

2. Continue to manipulate the lymph node groups one at a time and by moving towards the treatment site.

3. Treat the treatment site. (around 5 minutes)



5. End of treatment



End the treatment by manipulating the lymph node area of the venous angle.

4. Move back from the treatment site via the same lymph node groups and back to venous angle.

Small handles - great benefits!
6D TAPE®


Small handles - great benefits!

The lymphatic system – watershed areas and directions of flow




Small handles - great benefits!
6D TAPE®

Small handles - great benefits!



Six-Dimensional
Activation

The lymphatic system: venous angles




Most of the lymph is returned to blood circulation through the **venous angles (angulus venosus)**, which are junctions of about 1.5 to 2 centimeters in size at the bottom of the neck (in the middle of the supraclavicular fossae).

The left venous angle returns some three quarters of the body's lymph into circulation, and it covers the lymphatic fluid of the left upper quadrant of our bodies and both our lower limbs.

The right venous angle is responsible for returning about one quarter of the body's lymph back into circulation. It receives the lymph from the right upper quadrant of the body.

The lymphatic system – large lymph trunks

The body's largest lymphatic vessel is the thoracic duct (**ductus thoracicus**), which drains into the left venous angle.



6D Tape: Taping in Practice




6D Tape: Taping Instructions



1. Check – The skin must be clean, dry, and hairless.

2. Cut

- Cut the tape to the desired length and round the corners of the tape, or
- use a self-care kit with pre-cut tapes.



3. Tape

- Carefully remove the backing without touching the adhesive surface with your hands.
- The recommended stretch is 0%, but the tape can be stretched by 15% if necessary.
- Carefully rub the tape so that it adheres fully to the skin.

4. Activate – Use the handles to move the tissue in different directions: lift, twist, stretch and press.

5. Remove – The recommended wear time for 6D Tape is 2–10 days.

6D TAPE[®]
Six-Dimensional Activation

Small handles – great benefits!

Part 2: How to treat with 6D Tape?



Treatment protocol for 6DTape manipulation

Watershed areas and directions of flow

6D Taping in practice

Treatment videos / training with videos

summary

6D Tape is for active self-care, physiotherapy, and lymphatic therapy - 6D Tape is suitable for many scar, edema, and pain related indications



6D Tape® Self-Care Packages



6D Tape® 5-m Rolls



6D Action® Clinical Pro Negative Pressure Device Precisely Controlled Cupping

Six-Dimensional Activation
6D TAPE®

Small handles - great benefits!



6D TAPE[®] Six-Dimensional Activation
 Small handles – great benefits!

Active 6D
 tape with
 mobilization
 handles

For treating
**tennis
 elbow** and
golfer's elbow

Self-care package for tennis elbow and golfer's elbow. Includes instructions and precut 6D Tape strips for two 3-10-day treatments.


 Patented
 www.6dtape.com Copyright 6D Tape Inc. 2021

**Six-Dimensional
 Activation**

6D Tape self-care package

Tennis and Golfer's Elbow





Six-Dimensional
Activation

6D Tape – Universal tennis and golfers elbow




1. one minute
(10-15 x) left
neck

2. Manipulate elbow by starting lifting upwards from the first handles, then continue the middle handles and last the lowest handles.

3. While coming upwards use also direction pushing down and return the same way up.



6D Tape – Universal Tennis and Golfers Elbow



Six-Dimensional
Activation




6D Tape – Applications – never just one right solution



Recommended length to use is from 1 to 4 handles, in some special cases longer ones. Less restriction on joint or muscle function = feels more comfortable

The handle can be vertical or horizontal

Tape can be cut between the handles and the widthwise to a desired width, for example face and fingers is good $\frac{1}{4}$ of a size about 1 cm



6D Tape – Applications – never just one right solution



For lymphatic spider cut leave at least one handle uncut. The legs can be cut to 1-4

The lymphatic taping is not so restrictive in movements as uncut full handle therefore the lymphatic taping can be longer. The "legs" although can come off from the skin easier than full handle.

The smaller the treated area the smaller the 6D Tape should be.

Tape is not generally recommended to overlap but if so the ends should go at least one handle over to secure maximum adhesion



6D Tape – Wrist – Carpal Tunnel – De Quervain




The product packaging for 6D Tape is shown. It features the brand name "6D TAPE" in large blue letters, with "Six-Dimensional Activation" in smaller text above it. Below the name is the slogan "Small handles - great benefits!". A central image shows a hand with a piece of blue tape being applied. The text "Active 6D tape with mobilization handles" is prominently displayed. To the right, it says "For treating Carpal Tunnel and De Quervain Syndrome". Below this, smaller text reads: "Self-care package for carpal tunnel and de Quervain syndrome. Includes instructions and pre-cut 6D Tape strips for two 3-10 day treatments". At the bottom, there is a QR code, the word "Patented", the website "www.6dtape.com", the copyright notice "Copyright 6D Tape Inc. 2023", and the text "6D Tape self-care package".



Six-Dimensional Activation
6D TAPE[®]
Small handles - great benefits!

6D Tape – Wrist – Carpal Tunnel – De Quervain




Six-Dimensional
Activation




Knee treatments

Six-Dimensional Activation
6D TAPE[®]
Small handles – great benefits!

Six-Dimensional Activation
6D TAPE[®]
Small handles – great benefits!




Six-Dimensional Activation



Active **6D** for knee treatments
tape with mobilization handles

Self-care package for knee arthrosis, pain and swelling. Includes instructions and pre-cut 6D Tape strips for two 3-10 days knee treatments




Patented
www.6dtape.com Copyright 6D Tape Inc. 2021

6D Tape self-care package



Knee & Jumper's knee



6D Tape – Universal knee + Patellar Tendinitis



1. Manipulate one minute (10-15 x) left neck




2. Manipulate one minute (10-15 x) groin above the treated knee



3. Manipulate 5 minutes the treated knee.

4. Return the same way up. First groin (10-15 x) and last the left neck (10-15 x).



Patellar Tendinitis – jumpers knee

6D TAPE®
Six-Dimensional Activation

Small handles – great benefits!

6D Tape – Patellar Tendinitis – Jumpers' knee





6D Tape –Thigh muscles – Hamstrings



1. One minute in left venous angle, 2. one minute left groin, 3. Then 5 minutes on hamstrings and return one m. groin, one m. neck.

6D TAPE[®]
Six-Dimensional Activation

Small handles – great benefits!

Thigh – Quadriceps muscles

"By the end of the season I got 6D Tape for my use. This helped my recovery from an injury, alleviated pain, and relieved impingement caused by a scar," told Erik Back after finishing 400 meters with his personal best (47,14 s) in the contest between the national teams of Finland and Sweden.

National team athletes Erik Back and Eljas Aalto in the picture.




6D Tape – Calf muscles – Gastrocnemius & Soleus



1. Manipulate one minute (10-15 x) left neck




2. Manipulate (10-15 x) groin above the treated knee



3. Manipulate (10-15 x) the medial knee

4. Manipulate the calf from upper handles down and then up.

5. Return the same way up. (Knee, groin & neck)






6D Tape – Plantar Fasciitis – Heel Spur

The heel spur is recommended to 6D Tape with short 2 handle application which can be easily replaced compared to long taping. Secure the endings.

Handle can be cut into halves to secure better formation of the tape like here with achilles tendon.



Under the heel handle is cut off totally away so that it won't press inside the shoe.


6D Tape – Back of the foot



For the foot can be used single handles or three handles cut into the lymphatic pattern.



6D Tape – Upper back – Neck




Why 6D Tape is good for Wound and Scar Healing?



- **Enhanced Lymph Flow:** 6D Tape's unique ability to increase lymph flow is pivotal in accelerating wound healing and reducing edema.
- **Scar Tissue Management:** It effectively releases scar and fibrotic tissue, making it ideal for post-surgical recovery and C-Section scars.
- **Pain Relief:** The treatment provides significant pain relief by diminishing the swelling that presses on nerves, thereby enhancing patient comfort throughout the healing process

6D Tape – Shoulder



Six-Dimensional
Activation








6D Tape – Lower back – Hip



If one has stiff shoulder joint and difficult to reach on the back, the horizontal tape in the lower back can be from 5 to 7 handles.



6D Tape directly over the old scar



Usually 2-4 week after the operation you can go on top to the scar when it is totally healed and no possibilities of causing extra tissue damage.

Treatment is focused directly under the scar.

Especially small areas like fingers, toes, palm or face where space limitations for multiple strips.



6D Taping around the old scar

This approach offers a significant advantage as it enables multidirectional stretching of the scar tissue by applying force from different sides. Taping around the scar is particularly suitable under the following conditions:

1. Enhanced Stretching Force:

Leaving the scar uncovered while taping around it allows for stronger stretching forces to be applied directly to the scar tissue.

2. Adhesion Release:

Targeting tissues adjacent to the scar helps release adhesions, improving the mobility of connective tissues and enhancing overall flexibility.

6D Taping around the old scar



3. Pain or Discomfort with Direct Taping:

When direct taping and mobilization cause pain, irritation, or discomfort, taping around the scar provides a gentler alternative.

4. Alleviation of Referred Pain:

Releasing tension in the surrounding tissues can reduce discomfort or referred pain caused by scar-related adhesions and restricted connective tissue chains.



6D Taping both Directly over the scar and around the old scar


6D Tape – Post surgical scar therapy



Total thigh muscle rupture after slipping on a watery golf green. 9 weeks forbidden to bend the knee. Left pic. 3 days after operation. Very important to start mobilization on lymph nodes above the knee as soon as possible from the operation. (Venous angle and groin lymph nodes.) Recovery will happen much faster



6D Tape – Post surgical scar therapy



Six-Dimensional
Activation

Usually 2-4 week after the operation you can go on top to the scar when it is totally healed and no possibilities of causing extra tissue damage.



Post surgical scar therapy




Lymphatic taping
without going on top
of scar – too sensitive
still

Lymphatic taping
with 3 half handles
on top of scar – scar
not sensitive
anymore

Half handles
around the scar
and groin + inner
tight activation 2
handle strips

Combination
Lymphatic taping
with 6D Tape and
regular kinesiology t.



Our Customers Love Us



“We are using the manual version of 6D Tape in the UK. The feedback from the patients and practitioners has been fantastic.”

“Since the new 6D Tape® Pro automated tape will work 24/7, even the most severe cases can perhaps be managed without continuous bandaging or intensive therapy. This is another huge advantage to its use in our field”



- **Denise Hardy**, Clinical Nurse Specialist, Kendal Lymphology Centre, Kendal, UK

<https://www.kendal-lymphology.co.uk/>

Our Customers Love Us



“In our opinion the new automated 6D Tape® Pro therapy tape will be a revolutionary breakthrough in lymphedema, swelling and scar treatments.”

Patryk Gawrysiak, Specialist Physiotherapist in Lymphoedema/Team Leader, and PT Anca Tudor, St Georges Hospital, London, UK





Active 6D Tape with mobilization handles for Scars, Swelling and Pain in UK

➤ Premium tape 2 handles:

10 pcs x 5cm x 12cm, PIP code 434-4818

➤ Premium tape 4 handles:

5 pcs x 5cm x 24cm, PIP code 434-4800

➤ Premium tape:

2 pcs x 5cm x 1m, PIP code 434-4628

➤ Premium tape:

1 pcs x 5cm x 5m, PIP code 434-4602

6D Tape



Thank you for attending

Questions and comments



01/14/2025 8:55:53 AM

